

MULTIGATE CubeWare® 2.0



CUBEWARE® V'S CUBEWARE® 2.0 (NON-STERILE)

CubeWare®



PRODUCT DETAILS	Bowl (B) 400mL 108 x 108 x 72mm	Bowl (B) 1.0L 152 x 127 x 95mm	Bowl (B) 1.5L 175 x 150 x 103mm	Bowl (B) 2.4L 213 x 170 x 112mm	Bowl (B) 7.5L 366 x 284 x 129mm	Tray (B) 400mL 219 x 74 x 46mm	Tray (N) 700mL 253 x 89 x 52mm	Tray (G) 700mL 253 x 89 x 52mm	Tray (Y) 700mL 253 x 89 x 52mm	Tray (B) 1.0L 290 x 103 x 55mm	Tray (Y) 1.0L 290 x 103 x 55mm	Gallipot (C) 20-120mL 64 x 50 x 51mm	Gallipot (C) 50-250mL 91 x 50 x 62mm
PIECE WEIGHT	24.5 g	45.5 g	60.1 g	91.6 g	246 g	30.3 g	41 g	41 g	41 g	58 g	58 g	13.9 g	16.5 g
CODE	DP 1817	DP 1812	DP 1822	DP 1820	DP 1810	DP 1816	DP 1814	DP 1815	DP 1813	DP 1811	DP 1821	CW 1010	DP 1823
GHG RECYCLE*	70 gC02e	129 gC02e	174 gC02e	269 gC02e	729 gC02e	82 gC02e	115 gC02e	115 gC02e	115 gC02e	173 gC02e	167 gC02e	40 gC02e	46 gC02e
GHG LANDFILL*	106 gC02e	196 gC02e	263 gC02e	406 gC02e	1,107 gC02e	125 gC02e	176 gC02e	176 gC02e	176 gC02e	263 gC02e	253 gC02e	60 gC02e	71 gC02e
GHG INCINERATE*	168 gC02e	311 gC02e	416 gC02e	643 gC02e	1,756 gC02e	199 gC02e	280 gC02e	280 gC02e	280 gC02e	416 gC02e	400 gC02e	95 gC02e	113 gC02e

CubeWare® 2.0



PRODUCT DETAILS	Bowl (T) 400mL 131 x 109 x 60mm	Bowl (T) 1.0L 175 x 142 x 79mm	Bowl (T) 1.5L 199 x 165 x 89mm	Bowl (T) 2.5L 237 x 195 x 106mm	Bowl (T) 7.5L 360 x 277 x 150mm	Tray (T) 400mL 214 x 77 x 48mm	Tray (C) 700mL 255 x 94 x 55mm	Tray (T) 700mL 255 x 94 x 55mm	Tray (Y) 700mL 255 x 94 x 55mm	Tray (T) 1.0L 284 x 99 x 63mm	Tray (Y) 1.0L 284 x 99 x 63mm	Gallipot (C) 20-120mL 75 x 75 x 51mm	Gallipot (C) 50-250mL 96 x 96 x 63mm
PIECE WEIGHT	19 g	36 g	56 g	80 g	180 g	21 g	34 g	34 g	34 g	42 g	42 g	8 g	14 g
CODE	HAA-0017	HAA-0019	HAA-0020	HAA-0021	HAA-0025	HBA-0015	HBA-0017	HBA-0018	HBA-0019	HBA-0020	HBA-0021	HEA-0006	HEA-0007
GHG RECYCLE*	49 gC02e	95 gC02e	149 gC02e	207 gC02e	475 gC02e	54 gC02e	88 gC02e	88 gC02e	88 gC02e	108 gC02e	108 gC02e	20 gC02e	36 gC02e
GHG LANDFILL*	77 gC02e	148 gC02e	232 gC02e	325 gC02e	741 gC02e	85 gC02e	138 gC02e	138 gC02e	138 gC02e	170 gC02e	170 gC02e	32 gC02e	52 gC02e
GHG INCINERATE*	126 gC02e	239 gC02e	374 gC02e	528 gC02e	1,198 gC02e	138 gC02e	225 gC02e	224 gC02e	224 gC02e	277 gC02e	277 gC02e	52 gC02e	92 gC02e

TABLE KEY: (B) Blue, (C) Clear, (G) Green, (N) Natural, (T) Teal, (Y) Yellow, (GHG) Greenhouse Gas Emissions

Reduce environmental impact



Optimised to reduced its environmental impact — our latest innovation, CubeWare® 2.0 is an evolution in our hollowware line. With a focus firmly on the planet we have designed CubeWare® 2.0 to deliver impact reductions across all key indicators measured.



Less waste¹
CubeWare® 2.0 has been designed to achieve more with less. It uses up to 42% less plastic than our first generation CubeWare® hollowware.



Reduced impact²
Less plastic and a less energy intensive manufacturing process combine to significantly reduce CubeWare® 2.0's environmental footprint.



Recyclable³
CubeWare® 2.0 hollowware is made from 100% polypropylene (PP), a recyclable and common plastic used in hospital products.

PHYSICAL WASTE REDUCED BY UP TO:⁴

42%

21.8% AVERAGE WEIGHT REDUCTION⁷

CRADLE-TO-GRAVE GREENHOUSE GAS EMISSIONS REDUCED BY UP TO:⁵

46%

29.7% AVERAGE GHG EMISSIONS REDUCTION⁷

CRADLE-TO-GRAVE WATER USE REDUCED BY UP TO:⁶

55%

30% AVERAGE WATER USE REDUCTION⁷

For more information please contact your local Business Development Manager on (02) 8718 2888 or visit multigate.com.au