

Dressing Pack

08-888NP (Previous Version)



Basic Dressing Pack with Gauze Swabs and Cotton Balls

Re-order code: **08-888NP** (Previous Version)

Size	-
Units per Inner	20 pks
Inners per Carton	8 bags
Carton Quantity	160 pks
Sterile	Sterile
Use	Single use
Colour	-

Pack Contents

Description	Size	Quantity in Pack
Sterile Field	49cm x 49cm	1
Dressing Tray Blue	13.7cm x 13.7cm x 1.63cm	1
Paper Towel	40cm x 38cm	1
Non-Woven Swabs	7.5cm x 7.5cm	2
Cotton Balls		5
Plastic Forceps Yellow	10.8cm	1
Plastic Forceps Blue	10.8cm	2