

# All-Purpose Towels

## An economic alternative

A versatile towel that is highly absorbent with fast fluid capture rate. Can be used wet or dry.

## Features and benefits

High absorbency with quick saturation rate makes our All-Purpose Towels suitable for a wide range of cleaning applications.

- ✓ Absorbs up to 8.5 times its weight in fluid
- ✓ Fast fluid capture
- ✓ Superior wet and dry strength
- ✓ Low linting
- ✓ Latex free

### Multigate All-Purpose Towel - Product Specifications (4plys folded)

MATERIAL/ CONSTRUCTION	WEIGHT OF TOWEL	WATER ABSORBENCY	SINKING TIME	TENSILE STRENGTH: Machine Direc- tion	TENSILE STRENGTH: Cross Direction
50% Viscose 50% Polyester	45gsm	8.5 g/g	< 1.3 sec	98 N	28 N



### Multigate All-Purpose - Low Linting Towel

RE-ORDER	DESCRIPTION	SIZE (CM)	STERILE / NON-STERILE	PIECES PER DISPENSER	PIECES PER CARTON
57-885	All Purpose Towels Small	35 x 30	Non-Sterile	100pcs	1200pcs
57-886	All Purpose Towels Large	35 x 60	Non-Sterile	100pcs	1200pcs