

# MULTIGATE

## SafePro® 40

### SURGICAL GOWNS



## Protection and Comfort.

SafePro®40 surgical gowns offer superior protection and comfort. Their 4-layer polypropylene construction delivers an excellent fluid and microbial barrier. The lightweight comfortable fabric is perfect for instrument nurses, assistants and surgeons.

#### TECHNOLOGY

### Protection.

SafePro®40 gowns are manufactured with a high quality four-layer SMMS fabric that delivers a high-performance splash resistant barrier that repels fluids and bacteria.

#### FEATURES

### Comfort.

The SafePro®40 fabric is low linting, easily draped. Generous sleeve sizes and a range of lengths ensure the SafePro®40 range is designed with your comfort in mind.

### AAMI level 3 protection.

SafePro®40 standard and reinforced gowns meet and exceed AAMI level 3 standards making them the ideal solution for procedures involving light to moderate fluid levels.

### Reinforced.

SafePro®40 reinforced gowns feature a strong, yet breathable extra layer of reinforcement in the critical arm and front areas for added protection against fluid strikethrough.

- ✓ Comfortable
- ✓ Breathable
- ✓ Splash resistant barrier
- ✓ AAMI level 3 protection
- ✓ Reinforced critical areas
- ✓ Low linting



Microbial barrier that delivers improved barrier performance



Lightweight & breathable

#### Standard Gowns

RE-ORDER Gown +2 Towels	RE-ORDER Gown +1 Huck Towel	SIZE	GOWN LENGTH	STERILE / NON STERILE	PACKS PER CARTON	RE-ORDER Non-sterile Gown only
66-401	66-401H	Medium	112cm	Sterile	20 Packs	66-401NS
66-402	66-402H	Large	125cm	Sterile	20 Packs	66-402NS
66-403	66-403H	X-Large	135cm	Sterile	20 Packs	66-403NS
66-404	66-404H	XX-Large	135cm	Sterile	20 Packs	66-404NS

#### Reinforced Gowns

RE-ORDER Gown +2 Towels	RE-ORDER Gown +1 Huck Towel	SIZE	GOWN LENGTH	STERILE / NON STERILE	PACKS PER CARTON	RE-ORDER Non-sterile Gown only
66-411	66-411H	Medium	112cm	Sterile	20 Packs	66-411NS
66-412	66-412H	Large	125cm	Sterile	20 Packs	66-412NS
66-413	66-413H	X-Large	135cm	Sterile	20 Packs	66-413NS
66-414	66-414H	XX-Large	135cm	Sterile	20 Packs	66-414NS

For more information please contact your Business Development Manager:  
**(02) 8718 2888** or visit [multigate.com.au](http://multigate.com.au)