

MULTIGATE Softpro® Advance

SURGICAL GOWNS

Ultimate Protection and Performance.

Softpro® Advance gowns deliver ultimate protection and performance. The inner layer, made from an impervious membrane, provides compliance with AAMI level 4 standards, making them perfect for long duration, high fluid level procedures.



TECHNOLOGY

Breathable barrier.

Softpro® Advance gowns are made with a unique fabric that delivers a breathable viral barrier that is lightweight and comfortable under the most demanding operative conditions. They provide complete protection against blood and body fluid strikethrough and viral penetration.

FEATURES

AAMI level 4 protection.

Heat sealed seams and an impervious membrane ensures our Softpro® Advance gowns meet and exceed AAMI level 4 standards, making them the ideal solution for high fluid, long duration procedures.

Comfort.

With generous dimensions to support upper body movement and a lightweight back panel to assist in keeping you cool, Softpro® Advance gowns are designed to maximise comfort during procedures. The neckline also features a soft and adjustable closure to suit different neck sizes.

- ✓ Impervious
- ✓ AAMI level 4 protection
- ✓ Breathable
- ✓ Lightweight back panel
- ✓ Viral barrier
- ✓ Low linting

Standard Gowns					
RE-ORDER Gown +2 Towels	SIZE	GOWN LENGTH	STERILE/ NON-STERILE	PACKS PER CARTON	RE-ORDER Non-sterile Gown only
27-001SA	Medium	112cm	Sterile	20 Packs	27-101SA
27-005SA	Large	125cm	Sterile	20 Packs	27-105SA
27-006SA	X-Large	135cm	Sterile	20 Packs	27-106SA
27-007SA	XX-Large	135cm	Sterile	20 Packs	27-107SA
27-002SA	Large, X-Long	144cm	Sterile	20 Packs	27-102SA
27-003SA	X-Large, X-Long	146cm	Sterile	20 Packs	27-103SA
27-004SA	XX-Large, X-Long	146cm	Sterile	20 Packs	27-104SA



Lightweight and breathable



Protection against fluid strikethrough and viral penetration

For more information please contact your Business Development Manager:
(02) 8718 2888 or visit multigate.com.au