

# MULTIGATE

## SoftPro®

### SURGICAL GOWNS



## Advanced Protection.

SoftPro® surgical gowns offer superior protection and comfort. Their 5-layer polypropylene construction delivers an effective splash resistant barrier against bacterial and fluid transfer whilst maintaining breathability, ensuring you stay comfortable and cooler during OR procedures.

#### TECHNOLOGY

### Protection.

SoftPro® gowns are manufactured with a high quality five-layer SSMMS fabric that delivers a high-performance splash resistant barrier that repels fluids and bacteria.

#### FEATURES

### Comfort.

The SoftPro® fabric is low linting and easily draped. With generous sleeve sizes and range of lengths, the SoftPro® range has been designed with your comfort in mind.

### AAMI level 3 protection.

SoftPro® gowns meet and exceed AAMI level 3 standards making them the ideal solution for procedures involving light to moderate fluid levels.

### Reinforced.

SoftPro® reinforced gowns feature a strong, yet breathable extra layer of reinforcement in the critical arm and front areas for added protection against fluid strikethrough.

- ✓ Comfortable
- ✓ Breathable
- ✓ Splash resistant barrier
- ✓ AAMI level 3 protection
- ✓ Reinforced critical areas
- ✓ Low linting

Standard Gowns						
RE-ORDER Gown +2 Towels	RE-ORDER Gown +1 Huck Towel	SIZE	GOWN LENGTH	STERILE / NON-STERILE	PACKS PER CARTON	RE-ORDER Non-sterile Gown only
98-001	98-001H	Medium	112cm	Sterile	20 Packs	99-101
98-002	98-002H	Large	125cm	Sterile	20 Packs	99-102
98-003	98-003H	X-Large	135cm	Sterile	20 Packs	99-103
98-004	98-004H	XX-Large	135cm	Sterile	20 Packs	99-104
98-006	98-006H	XXX-Large	140cm	Sterile	20 Packs	99-106
98-005	98-005H	X-Large, X-Long	148.5cm	Sterile	20 Packs	99-105
98-022	66-413H	Large *2 gown pack	125cm	Sterile	10 Packs	99-102

\*Includes two large gowns and either 4 standard towels or 2 Huck Towels

Reinforced Gowns						
RE-ORDER Gown +2 Towels	RE-ORDER Gown +1 Huck Towel	SIZE	GOWN LENGTH	STERILE / NON-STERILE	PACKS PER CARTON	RE-ORDER Non-sterile Gown only
99-001	99-001H	Medium	112cm	Sterile	20 Packs	-
99-002	99-002H	Large	125cm	Sterile	20 Packs	99-202
99-003	99-003H	X-Large	135cm	Sterile	20 Packs	99-203
99-004	99-004H	XX-Large	135cm	Sterile	20 Packs	99-204
99-006	99-006H	XXX-Large	140cm	Sterile	20 Packs	99-206
99-005	99-005H	X-Large, X-Long	148.5cm	Sterile	20 Packs	99-205



Microbial barrier that delivers improved barrier performance



Lightweight & breathable

For more information please contact your Business Development Manager:  
**(02) 8718 2888** or visit [multigate.com.au](http://multigate.com.au)