

MULTIGATE

MediCel®



Absorbency and strength you can trust.

A versatile towel that is highly absorbent and can perform when wet or dry.

Features and benefits

Made from non-woven spun lace cellulose/polyester/viscose fabric, MediCel® is ideal for a wide range of cleaning applications. Its high strength, absorbency and lasting performance when wet makes it the perfect choice for high fluid procedures, such as mopping up spills and operating room hand drying.

- ✓ Absorbs up to 8 times its weight in fluid
- ✓ Superior wet and dry strength
- ✓ Last performance, even when wet
- ✓ Soft and gentle
- ✓ Low linting
- ✓ Latex free

MediCel® - Product Specifications (4plys folded)

MATERIAL/ CONSTRUCTION	WEIGHT OF TOWEL	WATER ABSORBENCY	SINKING TIME	TENSILE STRENGTH: Machine Direction	TENSILE STRENGTH: Cross Direction
50% Cellulose 35% Polyester 15% Viscose	56gsm	up to 8 g/g	< 5 sec	120 N	35 N



MediCel® - All Purpose Towel

RE-ORDER	DESCRIPTION	SIZE (CM)	STERILE / NON-STERILE	PIECES PER DISPENSER	PIECES PER CARTON
57-120	MediCel® Towel Small	35 x 30	Non-Sterile	100pcs	1200pcs
57-121	MediCel® Towel Large	35 x 60	Non-Sterile	100pcs	1200pcs

For more information please contact your local Business Development Manager:

(02) 8718 2888 or visit multigate.com.au